

FXD-900


Digital Radiology System for Veterinary



FXD-900 gives precision, instant radiographic images with crystal clarity. No more consumption and waiting required during x-ray exposure with FXD-900. FXD-900 will be a convenience for diagnostic capability in a veterinarian practice.

Also you could check the previews with a monitor instantly.





BLADE VET Software

- High resolution image
- DICOM 3.0 Compatible
- DICOM modality worklist SCU and SCP
- Tissue equalization
- Sharpening
- Auto-shutter

Specification

- **Detector**
 - *Detector Array* : Full field, Single VHD CCD
 - *Active Image Size* : 43cm x 43cm
 - *Number of Active Pixel* : 3,072X3,072 (9M)
 - *Spatial Resolution* : 3.5 lp/mm
 - **Tube** : Standard 1.0/2.0mm/140KHU
Optional 0.6/1.2mm/300KHU
 - **Manual Driven Four Leaves Collimator**
 - **Generator**
 - Rating output 32kW
 - Tube voltage kV 40-125 (1kV step)
 - mA Range : 10-400mA (17 steps)
 - mAs Range : 0.1mAs – 400mAs
- 

- Dimensions

